
PRÁCTICA V.2
SÍNTESIS DE UN DERIVADO DE SULFANILNAMIDA
SYNTHESIS OF A SULFANILAMIDE DERIVATIVE

Silvia Pérez Silanes

C.I.F.A. Universidad de Navarra. 31008 Pamplona
sperez@unav.es; cifa@unav.es

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EXERCISE V.2

SYNTHESIS OF A SULFANILAMIDE DERIVATIVE


Silvia Pérez Silanes

C.I.F.A. Universidad de Navarra, 31008 Pamplona, Spain

E-mail: sperez@unav.es; cifa@unav.es

INTRODUCTION: ANTIBACTERIAL AGENTS

The sulfonamides are a group of substances derived from *p*-aminobenzenosulfonamide or sulfanilamide, which was initially discovered as an active product formed in the organism after administration of an antibacterial agent, the dye protosil.


While the SO_2NH_2 group, in itself, is not essential for activity, it is necessary that the sulfur be directly linked to the benzene ring. The group NH_2 , in the para position, is an indispensable requirement for maintaining activity.

With regard to the mechanism of action, the sulfonamides are antibacterial compounds due to inhibition of the biosynthesis of folic acids and are diuretic agents due to inhibition of the carbonic anhydrase enzyme.

There are many sulfanilamide derivatives that are used in pharmacology. The first ones used in chemotherapy presented problems with respect to their elimination, derived from their lack of solubility; because of their pK_a (approximately 6 or lower in bacterial infections), a great number of them are found in a non-ionized form which, upon being insoluble, precipitate in water. In order to avoid this problem, the pK_a of the sulfonamides was reduced by introducing electron acceptors on the nitrogen of the sulfonamide group. Among these electron acceptors, the heterocyclic derivatives, such as sulfisoxazole, sulfadiazine, sulfamerazine, or sulfametazine, are less toxic.

In general, the sulfonamides are usually prepared by chlorosulfonation of benzene rings with chlorosulfonic acid, followed by aminolysis. The preparation of antibacterial sulfonamides, derived from 4-aminobenzenosulfonic acid, is similar, from the 4-acetamidobenzenesulfonyl chloride, which is obtained by reacting chlorosulfonic acid with acetanilide. From a synthesis point of view, the problem usually lies in the preparation of the amine, usually heterocyclic, which is used in aminolysis.

SYNTHESIS SCHEME**B: SULFONAMIDE****EXPERIMENTAL PROCEDURE**

The synthesis consists of two stages/phases.

1. Reaction of acetanilide with chlorosulfonic acid, giving rise to the formation of *p*-acetamidobenzenesulfonyl chloride (A).
2. Reaction of *p*-acetamidobenzenesulfonyl chloride (A) with diethylamine, under reflux, in order to give *N,N*-diethyl-*p*-acetamidobenzenesulfonamide, which after hydrolysis reaction, yields *N,N*-diethyl-*p*-aminobenzenesulfonamide (B).

1. PREPARATION OF *p*-ACETAMIDO BENZENESULFONYL CHLORIDE***EXPERIMENTAL PROCEDURE*****Before beginning, dry all glassware thoroughly**

Chlorosulfonic acid (0.15 mol) is added dropwise to dry acetanilide (0.03 mol). Stirring is maintained constantly during the addition in order to homogenize the mixture.

Once the addition is completed, the mixture is heated to reflux for 1 h in order to complete the reaction. The mixture is allowed to cool and is then poured over crushed ice. The mixture is then stirred continuously until a suspension of a white solid is

obtained; this solid corresponds to *p*-acetamidobenzenesulfonyl chloride. This is filtered under vacuum, and the filtrate is washed with cold water. The filtrate is then dried thoroughly.

QUESTIONS

1. Why is it necessary that the material be well dried?
2. Why is the reaction mixture poured over ice?
3. What is the mechanism of reaction of sulfonation?

2. PREPARATION OF *N,N*-diethylbenzenesulfonamide

EXPERIMENTAL PROCEDURE

p-Acetamidobenzenesulfonyl chloride is placed in a flask and a mixture of 0.06 moles of diethylamine in 30 mL of H₂O is added. Mix well and stir while heating to just below boiling point, until the corresponding *N,N*-diethyl-*p*-acetamidobenzenesulfonamide (reaction observed using TLC) is obtained.

Cool the suspension over ice and then add diluted sulfuric acid until the mixture reaches an acidic pH. Filter the product under vacuum, wash the filtrate with cold water, and then dry completely.

For hydrolysis, heat *N,N*-diethyl-*p*-acetamidobenzenesulfonamide with 20 mL of water and 10 mL of concentrated hydrochloric acid (the reaction is observed using TLC). Once the reaction has ended, the content of the flask is poured into a beaker. A saturated solution of sodium bicarbonate is carefully added in small portions and with stirring. Check the pH of the suspension. If it is still acidic, add more sodium bicarbonate, until a neutral pH is reached. Cool with ice, filter the sulfonamide under vacuum, and dry. Recrystallize from the appropriate solvent.

INSTRUCTIONS AND QUESTIONS

1. Prepare the material and calculate the necessary quantities of commercial reagents.
2. The reactions are observed using TLC with a MP (mobile phase) of ethyl acetate/petroleum ether (1:1).
3. Once the product is obtained, the following data are noted:

Time of reaction of synthesis of *N,N*-diethyl-*p*-acetamidobenzenesulfonamide.

Time of reaction of hydrolysis.

M.P. (melting point)

IR-reference (acetanilide)

IR-product obtained (sulfonamide)

NMR (DMSO-d₆; 200 MHz) (δ):

ELEMENTAL ANALYSIS (C, H, N)


4. Briefly answer the following questions:

- What happens to the chlorine atom from the *p*-acetamidobenzenesulfonyl chloride?
- Why do we add sulfuric acid?
- Why do we add sodium bicarbonate at the end of the hydrolysis?
- What is the appropriate recrystallization solvent/solvent for recrystallization?

PREPARATION OF A REISSELT DERIVATIVE FROM 6-NITROQUINOLINE

INTRODUCTION: REISSELT REACTION

The chemistry of the Reissert compounds has been the object of numerous revisions since 1905. A. Reissert made the first example known: the reaction between quinoline (1 mol) and benzoyl chloride (2 mol) in an aqueous potassium cyanide medium.

CLASSIC EXAMPLE OF REISSELT
QUINOLINE DERIVATIVE

Numerous modifications of the reaction, in the substrate (generally, quinoline and isoquinoline) as well as in the cyanide-donating/cyanide-releasing agents and the reaction medium. Therefore, originally, the Reissert derivatives were prepared in an aqueous medium by slowly adding acid chloride (2 mol) to a suspension of quinoline (1 mol) in aqueous potassium cyanide (3 mol).

This method has numerous advantages:

- Neither the quinoline nor the reaction product is soluble in water.
- When the reaction is conducted in aqueous medium, a *pseudobase* is frequently isolated by adding an OH group instead of a CN group.
- In these conditions, the acid chloride is hydrolyzed.


The introduction of nonaqueous substances facilitated dilution and/or dissolution of the acid chloride as well as the substrate.

At that time, the reaction was carried out by slowly adding acid chloride (frequently dissolved/diluted in dry dichloromethane) to the nitrogenated base and the potassium cyanide in a dichloromethane medium to which the minimum quantity of water necessary for solubilizing the salt had been added.

Subsequently, better reaction yields were obtained using phase transfer agents, such as benzyltrialkylammonium chloride. With these agents, a decrease in the formation of the pseudobase was also obtained.

One of the most efficient cyanide-donating/cyanide-releasing agents, which permits the use of an anhydrous medium, is **trimethylsilane cyanide**. In this case, aluminum trichloride (a Lewis acid) is used as the catalyst.

A probable mechanism for the reaction (proposed by McEwen and Cobb) involves the formation of an intermediate: benzoylquinoline chloride. Subsequent addition of the cyano group to this intermediate is verified with preference in position 2 (see scheme below):

REISSERT REACTION MECHANISM


INESTABLE = UNSTABLE

The absence of water in the medium prevents hydrolysis of the more reactive acid chlorides, a competitive process which led to negative results in some of the attempts to carry out a reaction in an aqueous medium or in combination of the two phases, aqueous and organic.

SCHEME OF SYNTHESIS

EXPERIMENTAL PROCEDURE

A solution of acid chloride⁽¹⁾ (1.5 mmol) in dry dichloromethane (20 mL) is slowly added to a mixture of 6-nitroquinoline (1 g), trimethylsilane cyanide (1.5 mmol), and aluminium chloride (catalytic quantities) in dry dichloromethane (50 mL), under nitrogen atmosphere, and in an ice bath.

REISSERT DERIVATIVE OF
6-NITROQUINOLINE

After being observed by TLC and stirred at room temperature, the reaction mixture is poured onto water (50 mL) and washed with water (2 x 50 mL), HCl (dil.) (until the aqueous phase remains colorless), water (2 x 50 mL), aqueous solution saturated with sodium bicarbonate (until the aqueous phase remains colorless), and water (2 x 50 mL). The organic phase is dried over anhydrous sodium sulfate and filtered. The solvents are eliminated under reduced pressure.

The residue is purified by column chromatography. The solvents are removed under reduced pressure.

INSTRUCTIONS AND QUESTIONS

- 1. Prepare the material and calculate the necessary quantities of the commercial reagents.**
- 2. Reaction. The reaction is observed by TLC with a MP toluene/ethyl acetate (50:50) against the starting product, 6-nitroquinoline.** The reaction time in the given the work conditions will be determined.
- 3. Characterization of 6-nitroquinoline (starting product).**

IR-reference (6-nitroquinoline)

NMR H¹ spectra of 6-nitroquinoline is studied and each signal is assigned to the corresponding H. (See appendix.)

- 4. Once the product has been prepared, the following data are noted:**

Time of reaction

M.P. (melting point)

IR (Reissert derivative obtained)

(1): The lab monitor will indicate the work reagent to each group.

NMR (DMSO-d₆; 200 MHz) (δ) (Reissert derivative obtained):


ELEMENTAL ANALYSIS (C, H, N)

5. Briefly answer the following questions:


- What is the reaction product of trimethylsilane cyanide with water?
- How is the residual trimethylsilane cyanide eliminated?
- Could ethanol be used as a reaction solvent instead of dichloromethane?

APPENDIX: CHARACTERIZATION OF 6-NITROQUINOLINE BY RMN-H¹ SPECTROSCOPY


6. Nitroquinoline (chemical structure)

6-NITROQUINOLINA, ASPECTO GENERAL

Double resonance spectra for determining coupling/couplings:

Silvia Pérez Silanes
sperez@unav.es; cifa@unav.es

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, solvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

A cyanide antidote kit (amyl nitrite, sodium nitrite and sodium thiosulfate) should be available in any cyanide work area before work commences.

Chlorosulfonic acid is corrosive, reacts violently with water and is poisonous.

This document has been supervised by Prof. Silvia Pérez Silanes (sperez@unav.es) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.