PRÁCTICA V.9
PRACTICA DE ANTIHISTAMÍNICOS
ANTIHISTAMINES

Dra. María Elena Montoya Alfaro

Profesora Principal, T.P. (Cod. 04573A)

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ANTI-HISTAMINES


María Elena Montoya Alfaro and Eriberto Carrasco

*Faculty of Pharmacy and Biochemistry, Academic-Professional School of Pharmacy and
Biochemistry, San Marcos Major National University*

E-mail: MEMONTOYAA@latinmail.com

OBJECTIVE

To identify and assay H1 antihistamine drugs

CHLORPHENIRAMINE

MATERIALS

- Test tubes
- Mortar
- Dropping funnel
- Chromatographic microcamera
- Chromatographic plates
- Flat-bottom flasks, burette

REAGENTS

- Solvents: alcohol, ether, and chloroform
- Dragendorff reagent
- Iodine fumes
- Mayer reagent
- Iodine reagent
- Picric acid reagent
- Crystal violet TS
- Acetic acid
- Acetic anhydride
- Perchloric acid 0.1 N

SAMPLES

- Pure drug of reference or standard: chlorpheniramine maleate
- Pharmaceutical form or problem sample: 4-mg tablets of chlorpheniramine maleate

PROCEDURE

Preparing the problem sample:

- In a mortar, grind a tablet until it is a fine powder.
- Put powder in a test tube, add 3 mL of distilled water, shake strongly, and filtrate. Save the liquid phase for assays.

Identification:

Aspect, color, and smell of reference sample.

Solubility: Prove pure drug solubility in water, alcohol, chloroform, and ether

Thin-layer chromatography:

- For saturating the chromatographic tank, before 30 min, cover the walls with filter paper and add (mobile phase) solvent approximately in 1 cm of height.
- Prepare the sample and reference solutions approximately 1 % in the solvent of greater solubility and apply with a capillary tube at 0.5 cm of inferior and lateral margins. Apply sample, standard, and a sample/standard mixture.
- Let it dry and mark with a pencil the superior margin of solvent front (mobile phase). Put the chromatographic plate vertically inside the chromatographic tank and rapidly put on the lid.
- Let the solvent reach the marked superior border. Take the plate out of the tank and let it dry.
- Detect the stains by UV wavelength and/or reveal with Dragendorff or iodine reagents and encircle the stains.
- Determine the center of the stain and measure the covered distance by sample and standard, solvent front, then calculate R_f of the sample and the reference or standard and compare.

The R_f of a compound is a physical constant, that can vary with any physicochemical change, thus, it should always compare with a chromatographic standard in the same experimental conditions.

Identification reactions:

Mostly, they are based on the presence of the amino group in the side chain.

Do the assay with the problem sample and standard at the same time by comparing results.

- With picric reagent: To 1 mL of sample, add drops of reagent. You should observe a yellow precipitate.
- With iodine reagent: To 1 mL of sample, add drops of reagent. You should observe a brown precipitate.
- With Mayer reagent: To 1 mL of sample, add drops of reagent. You should observe a white precipitate.

Chlorpheniramine maleate assay:

It uses the method of base titration in anhydrous medium.

Weigh an amount of approximately 40 mg of the substance (the work group will receive the sample for evaluating already weighed, but they will not know the exact weight, which is what they will have to determine). Add 10 mL of acetic acid, 3 mL of acetic anhydride, and drops of crystal violet TS as indicator. Assay with 0.1 N perchloric acid. The end of the reaction is determined by change of color from blue to greenish blue. Do a blank of reagent if necessary. 1 mL of 0.1 N perchloric acid equals 19.54 mg of chlorpheniramine maleate.

PRACTICAL REPORT

TABLE NUMBER:

Date:.....

Table members:

- 1.
- 2.
- 3.
- 4.

Practice number:

Title:

Objectives:


RESULTS

IDENTIFICATION

Characteristics: 1.
 2.
 3.

Solubility: with solvent: A. B. C.

Scheme of obtained chromatogram:

 R_f : Standard: R_f : Problem sample

ASSAY

Obtained waste:

Calculations:

CONCLUSIONS:

OBSERVATIONS/RECOMMENDATIONS (from student):

María Elena Montoya Alfaro

MEMONTOYAA@latinmail.com

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

Picric acid reagent must be kept wet at all times since it may explode if allowed to dry.

This document has been supervised by Profs. María Elena Montoya Alfaro (memontoyaa@latinmail.com) and Eriberto Carrasco Raymunde (ercara55@yahoo.es) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.