PRÁCTICA VI.4

ANALISIS CUANTITATIVO DE MENTOL EN ACEITE ESENCIAL DE MENTA QUANTITATIVE ANALYSIS OF MENTHOL IN MENTHA ESSENTIAL OIL BY GAS CHROMATOGRAPHY

Morayma Molina, Isabel Cabello, Olga Lock*

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EXERCISE VI.4

QUANTITATIVE ANALYSIS OF MENTHOL IN *MENTHA* ESSENTIAL OIL BY GAS CHROMATOGRAPHY

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
**E-mail: olock@pucp.edu.pe*

1. INTRODUCTION

"Mint oil" is obtained from several mint species: *Mentha arvensis*, *M. spicata* (sin. *M. viridis*), *M. piperita*, *M. pulegium*, among others. Nevertheless, its main constituent is different in all of them, for example, *M. arvensis* oil contains 70–90 % of menthol, whereas in *M. spicata* predominates (–)-carvone and in *M. pulegium*, pulegone; this is explained because of the easy monoterpenes interconversion.

Menthol is a monoterpenes as crystalline needles having the following properties: m.p. 41–43 °C; b.p. 212 °C; d 0.89 g/mL (10 % alcoholic solution). It is slightly soluble in water and very soluble in alcohol, chloroform, ether, and petroleum ether.

Both the mint oil and menthol are extensively used in pharmaceutical preparations, mainly in the treatment of respiratory problems. Other uses are as a flavor ingredient and fragrant component. Menthol is also used as carminative and gastric sedative.

This experiment will describe the quantitative analysis of menthol in mint oil using the gas chromatography technique.

EXPERIMENTAL PROCEDURE

1.1. Essential oil extraction

- Weigh 100 g of dried crushed mint leaves and obtain the essential oil by steam distillation.

- b. Separate the distilled essential oil and add anhydrous sodium sulfate. Let it dry for several minutes, then remove the oil and estimate the yield
- c. Obtain the IR spectrum and analyze the main absorptions.

1.2. Gas chromatography (GC) analysis

1.2.1.External standard method

- a. Place 10 mg (or enough quantity in order to obtain the concentration of the menthol within the range of the calibration curve) of the essential oil in a 10-mL volumetric flask and bring up to volume with ethanol (M solution).
- b. Inject 1 μ L in the GC equipment, under the following conditions:
 - Column: HP-5, 30 m x 0.32 mm x 0.25 μ m (5 % phenylmethylpolysiloxane)
 - Injector temperature: 200 °C
 - Detector temperature: 300 °C
 - Oven temperature: 110 °C
- c. Prepare the calibration curve with ethanolic solutions containing 1.0, 3.0, 5.0, 7.0, and 10.0 mg/mL of menthol standard.
- d. Calculate using the calibration curve, the % of menthol, m, in the M solution, expressed in mg/mL.
- e. Calculate the % of menthol in the mint leaves, with the formula:

$$\% \text{ menthol} = 10 \times m/W \times \% \text{ st}$$

m: mg of menthol in M solution

W: mg of sample (dried mint leaves)

% st: purity of standard menthol

1.2.2.Internal standard method

- a. Weigh, separately, 10 mg of standard menthol and 10 mg of octanol (internal standard) in 10-mL volumetric flasks and bring up both to volume with ethanol (A and B solutions, respectively).
- b. Place 1 mL of A solution and 1 mL of B solution, in a 5-mL volumetric flask, bring up to volume with ethanol (C solution).

- c. Weigh 12–15 mg of mint oil in a 10-mL volumetric flask, bring up to volume with ethanol (D solution).
- d. Place 1 mL of D solution and 1 mL of B solution in a 5-mL volumetric flask, bring up to volume with ethanol (E solution).
- e. Inject 1 μ L of solutions C and E separately in the gas chromatograph, as described in 2.2.1 method conditions.
- f. Calculate the % of menthol in the mint oil with the formula:

$$\% \text{ of menthol} = 5 \times (A_m/A_{\text{ist E}}) / (A_{\text{st}}/A_{\text{ist C}}) \times W_{\text{st}}/W_m \times \% \text{ st}$$

A_m : menthol area in the sample, E solution

$A_{\text{ist E}}$: octanol (internal standard) area, E solution

A_{st} : menthol standard area, C solution

$A_{\text{ist C}}$: octanol (internal standard) area, C solution

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Olga Lock

olock@pucp.edu.pe

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories.

The handling of electrical instruments, heating elements, glass materials, solvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. OLGA LOCK SING (olock@pucp.edu.pe) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.