

**DIVISION OF ANALYTICAL CHEMISTRY  
COMMISSION ON MICROCHEMICAL TECHNIQUES\***

**RECOMMENDED TEST SUBSTANCES FOR THE  
MICRODETERMINATION OF OXYGEN IN  
ORGANIC COMPOUNDS**

A number of compounds are recommended for use as reference substances for the microdetermination of oxygen in organic substances. All of these substances, or a proper selection from the list, may be used to determine the universal applicability of a given method, either already described or one which might be developed in the future.

The compounds selected are stable over long periods of time and are non-hygroscopic (any exceptions to the latter generalization are so noted). The substances are either commercially available in a sufficiently pure state to be used for test purposes based on the accuracy of present-day methods, or may be purified or prepared by conventional laboratory means to meet these standards.

The compounds selected include the following:

- (1) compounds for the determination of the blank value and those which have extreme values (high as well as low) for oxygen content;
- (2) compounds representing a variety of structural types;
- (3) compounds containing elements which may cause interferences in the oxygen determination.

These recommendations may be changed or supplemented as the need arises.

\* Al Steyermark, Chairman (U.S.A.), W. Schöniger, Secretary (Switzerland), A. A. Benedetti-Pichler (U.S.A.), N. E. Gel'man (U.S.S.R.), S. A. Hölzel (Austria), G. Ingram (U.K.), J. Körbl (Czechoslovakia), R. Lévy (France), A. J. Llacer (Argentina), R. J. Magee (U.K.), H. Spitzky (Austria), M. Věčeřa (Czechoslovakia).

## SECTION OF ANALYTICAL CHEMISTRY

<i>List of substances</i>	<i>Empirical formula</i>	<i>Molecular weight</i>	O (%)	C (%)	H (%)	N (%)	Hal (%)	S (%)	<i>Element (%)</i>
Blank (C, H, N)									
Anthracene	$C_{14}H_{10}$	178.234		94.34	5.66				
Naphthalene	$C_{10}H_8$	128.174		93.71	6.29				
Diphenylamine	$C_{12}H_{11}N$	169.227		85.17	6.55	8.28			
$C, H(O)$									
Benzoic acid	$C_7H_6O_2$	122.123	26.20	68.85	4.95				
Cholesterol	$C_{27}H_{46}O$	386.664	4.14	83.87	11.99				
d-Glucose (dextrose)	$C_6H_{12}O_6$	180.156	53.28	40.00	6.71				
$N(C, H, O)$									
Acetanilide	$C_8H_9ON$	135.166	11.84	71.09	6.71	10.36			
2,4-Dinitrophenylhydrazine	$C_6H_6O_4N_4$	198.138	32.30	36.37	3.05	28.28			
$F(C, H, O \dots)$									
<i>p</i> -Fluorobenzoic acid	$C_7H_5O_2F$	140.113	22.84	60.01	3.60		13.56 F		
Trifluoroacetanilide	$C_8H_6ONF_3$	189.136	8.46	50.80	3.20	7.41	30.13 F		
$Br(C, H, O \dots)$									
<i>p</i> -Bromoacetanilide	$C_8H_8ONBr$	214.067	7.47	44.89	3.77	6.54	37.33 Br		
$Cl(C, H, O \dots)$									
Chloroacetamide	$C_2H_4ONCl$	93.513	17.11	25.69	4.31	14.98	37.91 Cl		
$I(C, H, O \dots)$									
<i>o</i> -Iodobenzoic acid	$C_7H_5O_2I$	248.019	12.90	33.90	2.03		51.17 I		

